中国科学院水生生物研究所  
INSTITUTE OF HYDROBIOLOGY, CHINESE ACADEMY OF SCIENCES

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
The screenshot displays a software interface for X-ray data reconstruction, featuring several overlapping windows:

- Reconstruction (Top):** A window with tabs for Start, Settings, Advanced, Output, and Summary. It contains a plot showing a sharp peak on a baseline.
- Reconstruction (Right):** A smaller window showing the 'Advanced' settings tab, including a 'Position' field set to 1263, a 'fastest' field set to 1263, and a 'Top' field set to 2591.
- Reconstruction (Bottom Right):** A window showing the 'Advanced' settings tab with various options checked:
  - Smoothing (3)
  - Misalignment compensation:
    - Show -11.5
  - Object larger than field of view
  - Ring artifacts reduction (22)
  - Beam-hardening correction (30%) (prev: 30%)
  - CS rotation (deg) 0.00

- Batch manager (Bottom Left):** A window with buttons for Start batch, Submit again, Move up, Remove, and Remove all, along with a table for Batch job, Status, and Message.
- Color (Bottom Right):** A window with an 'Inverse' checkbox and a color gradient bar.
- Background:** A dark image showing a circular object, possibly a slice of a specimen, with a green circle overlaid.

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50  
455

```
image crop origin X=0  
Image crop origin Y=0  
Camera binning=2x2  
Image Rotation=0.2390  
Gantry direction=CW  
Number of connected scans=1  
Image Pixel Size (um)= 17.76  
Object to Source (mm)=120.625  
Camera to Source (mm)=171.118  
System Matrix Calibration=NO  
Vertical Object Position (mm)=60.869  
Optical Axis (line)= 653  
Filter=Al 0.5mm  
Image Format=TIFF  
Data Offset (bytes)= 264  
Horizontal overlap (pixel)=64
```

CTion

Ionizing xray dose calculator

Xray filter: 0.5mm Al

Applied voltage: 50  
Available range: 40 - 80 kV

Source current (μA): 455

Source-object distance (mm): 120.625

Dose Rate (mGy/min)

**BRUKER**

Air: 333.09  
Mouse (25mm diam.): 153.37  
Rat (75mm diam.): 57.22

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Calculate Save

Type Of Motion=STEP AND SHOOT  
Source Temperature=24.80 °C  
Study Date and Time=Oct 10, 2019 10:30:54  
Scan duration=00:04:45  
[Reconstruction]  
Reconstruction Program=NRecon

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